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Crystallization is used to purify a solid. The process requires a suitable solvent. A suitable solvent is one which readily dissolves the solid (solute) when the solvent is hot but not when it is cold. The best solvents exhibit a large difference in solubility over a reasonable range of temperatures. (eg, Water can be a crystallization solvent between 0-100°C; hydrocarbon solvents such as hexanes or petroleum ether have a different T range since they can be cooled below 0 degrees but boil below 100 degrees).

Characteristics of a solvent:

- chosen for solubilizing power-- solubility usually increases with increasing T
- polarity is important--like dissolves like; polar compounds are more soluble in polar solvents; nonpolar compounds in nonpolar solvents
- should be INERT but few are; eg, acetic acid is sometimes used as a solvent although it will certainly react with basic compounds
- almost all solvents are COMBUSTABLE--avoid flames
- mixed solvents (eg; 1:1 water/methanol) provide a huge range of possible solvents but they must be soluble in one another

Q: Is 95% ethanol a mixed solvent?

20/12/21, 09:36

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 2

Use solvent to get solids into solution but to get them out of solution:

- lower the temperature--solute will be less soluble
- concentrate the solution by removing solvent with a hot plate, heating mantle (flasks), steam bath (use in hood) or with the Roto-Evaporator.

To remove solvent:

1. You must have ebullition to concentrate at atmospheric pressure--use a boiling stone, a capillary tube, or agitation.
2. If you used reduced pressure to concentrate solution, use the water aspirator with a TRAP in the line. DO NOT turn off the water until the pressure is released. In general, CLAMP any flask that could conceivably trip over.
3. Do NOT use ebullition if using the Roto-Evaporator. The rotation provides sufficient agitation.

20/12/21, 09:37

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 3

Recrystallization or crystallization

- use an Erlenmeyer flask, it is specifically designed for this purpose
- dissolve solid in minimum amount of boiling solvent - add solvent in small amounts. For example, if you add 5 mL and approx. half of the solid dissolves, it should take only another 5 mL to dissolve the remaining half. If some of the solid does not dissolve then....
- are remaining particles your compound or insoluble material (eg, sand, old boiling stones)?
- to determine this, add ca. 10% more hot solvent. If insoluble material, you can decant (carefully transfer solution into another flask leaving the insoluble material behind) or filter.
- if filtering is necessary, do so to remove suspended solids, the faster the better, keep solution warm so crystallization does not occur (this may require filtering on a hot plate or other heating device).
- to "decolorize", use a small amount of charcoal and filter with "filter aid" (see below). For both (e) and (f), rinse filter paper with a small amount of hot solvent.
- let the filtered liquid (filtrate) cool to room temperature slowly in the Erlenmeyer flask
- cool the filtrate in an ice-water bath
- if crystals have not formed

DECOLORIZATION

Decolorization: Most organic compounds are colorless. Highly conjugated compounds (eg, polar polymers) will absorb light in the visible region of the spectrum and thus be "colored". If these highly polar, large molecules are impurities, they can be removed by use of finely granulated activated charcoal (Norit). Polar compounds (eg, polar impurities) adsorb to the charcoal which is insoluble in the solvent and can be filtered away from solution. Unfortunately, some of your compound will also adsorb if there is enough charcoal so the trick is to use just the right amount. Usually, a very small amount of charcoal will suffice (there is a lot of surface to these particles). The Norit is added in small amounts to the hot (but not boiling) solution until sufficient decolorization has occurred.

20/12/21, 09:38

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 5

CAUTION--trapped air in the Norit can cause rapid frothing when it hits the hot solution. The Norit can be filtered from the hot solution using fine filter paper or a **filter aid (Celite)** which is spread on top of the filter paper. To do this, make a slurry of the Celite in any solvent. Wet the filter paper and apply suction to make it stick. Now carefully pour the Celite slurry onto the filter paper so it evenly covers it. Then apply the suction. You may discard the solvent from the filter flask or change filter flasks prior to filtering the Norit containing solution.

20/12/21, 09:39

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 6

FILTERING

Used to remove insoluble solids suspended in solution.

Use a **GRAVITY FILTER FUNNEL** when you DON'T want the solid

Use a **HIRSCH** or **BUCHNER** funnel with vacuum when you do want the solid but...

Never use a Hirsch or Buchner funnel with a hot solution unless suggested by your instructor.

For the Hirsch and Buchner funnels, use a piece of round filter paper which fits the funnel. A proper fit is a piece of paper that just covers the holes but does not touch the sides of the funnel (ask for a demo). You will need to use reduced pressure, usually via a water aspirator (be certain that the TRAP is clean). For aqueous solutions, you can use the house vacuum line (be certain that the TRAP is clean).

20/12/21, 09:39

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 7

HEATING

There are different methods used for heating material in the laboratory. Flames are never used in the laboratory except in controlled situations (eg, isolated in fume hoods). Electric hot plates and heating mantles are most commonly used. Be careful not to turn this equipment to its highest setting which can burn it out. It does take several minutes for these instruments to reach the desired temperature. The heating mantles are plugged into a variable rheostat which provides a temperature control. Heating mantles are used for round-bottom flasks (rbf); choose an appropriate size to fit the flask you plan to use.

Steam is often used for heating volatile, non-aqueous, flammable solvents in°C maximum temperature from steam is sufficient for most commonly used organic solvents including many mixed aqueous solvents.

20/12/21, 09:40

COMMON LABORATORY TECHNIQUES - Chemistry LibreTexts

Serie 8

MELTING POINT

Melting Point Determination:

The standard physical property of a solid is its melting point. The melting point is actually a melting point range. It is used to help determine the purity of a solid and to help verify the identity of the compound. A pure compound should melt over a narrow temperature range. Impurities usually cause the melting point range to widen and lower in value. To obtain the melting point range, you record the temperature at which the first crystals begin to melt (solid to liquid phase) and the temperature at which the last crystal melts, eg; m.p. 124-126°C. Important considerations include: 1. do you have a representative sample of your compound or have you fished out the best looking crystal for a melting point; and 2. is the thermometer you are using correctly calibrated. If not, you will obtain incorrect values.

Extraction

Extraction is a method for moving a compound from one medium to another. For example, if you make coffee from coffee beans, you are extracting some flavorful components of the bean and some caffeine into the water. The remainder of the beans (grounds) are left behind and discarded. This is called a solid-liquid extraction. If you are trying to move a compound from one liquid phase (solvent 1) into another liquid phase (solvent 2), this is liquid-liquid extraction but the two solvents must be immiscible or insoluble to the extent that they form two distinct layers (why?). The compound is now distributed into two solvent layers which can be separated. By measuring the concentrations of compound in the two solvents (c_1 and c_2) we obtain a distribution coefficient, K , which is a constant for a given compound and given solvents *at a given temperature*, irrespective of the amounts of solvent present. A simple equation shows the relationship $K=c_2/c_1$